

Artificial Vision


Causes of blindness in the adult

- The most common causes of blindness
 - Age-related Macular Degeneration, AMD(나이관련황반변성)
 - Diabetic Retinopathy(당뇨망막병증)

- The most serious cause of blindness
 - Retinitis Pigmentosa (망막색소변성)

Care for Blindness


- visual impairment :
 - Legally considered percentage of lost body function: 24%(**Unioocular**) and 100%(**binocular**)
- 30% of blindness in the adult
 - Retinitis Pigmentosa **망막생소변성**; 1/4000(Normal)
 - AMD **나이관련황반변성**; 1/20(>aged 65)
- Artificial Organs in the eye :
 - Corneal clouding **각막혼탁** → Artificial Cornea
 - Lenticular opacity **수정체 혼탁** → Artificial Lens
 - Retina Damage → Artificial Retina?



BME2009



Structure of the eye

- retina(**망막**): nervous tissue
- choroid(**맥락막**): vascular tissue
- sclera(**공막**): connective tissue
- optic nerve(**시신경**)



Retina


- Recognizes image and transmits it to brain
- Can be compared to camera film



BME2009

Structure of the Retina



- photoreceptors(시세포층): Light Signal -> Electrical Signal
- RPE(망막색소상피층): Light absorption, Supply of Nutrition, a protective barrier
- Bruch's Membrane(브루크막) – Boundary tissue between retina and choroid



BME2009


Macula (황반, 黄斑)

- macula : center of retina, 5mm diameter, provides central vision
- Tinted yellow due to rich concentration of Xanthophyll(엽황소)



BME2009


The world that patients with macular degeneration see



BME2009

Recovery from retina damage

- Drug treatment > Target Cell issue
- Stem Cell > Differentiation issue
- Transplantation of visual cells
> Settlement issue
- Mimicking visual stimulation
with electrical stimulation




Why Artificial Retina?

- Progress for the last 20 years
 - Good preservation of the retina inner layers in the RP patients have become known.
 - success of cochlear implant
 - development in semiconductor technology
→ fabrication of a large-scale integrated circuit
 - AR research began in late 1980s

What happens when retina is electrically stimulated?

- Foerster (1929): patient saw a small spotlight on electrical stimulation of visual cortex - 'phosphene(인광)'
- Brindley & Lewin(1968)
 - Implant a visual cortex stimulator in a 52 years old patient
 - couldn't find conditions to obtain continuous images
- Dawson & Radtke (1977): electrical stimulation of retina
→ discovery of sense of phosphene



BME2009

Artificial vision

- Basic Concept:
replace function of retina cells with an electrical device.
- Hypothesis :
electrical stimulation can induce vision.
- Goal:
implantation of microchip for evoking vision.


BME2009

How many electrodes are required?

- Cochlear Implant
 - 6 electrodes \ll 30,000 auditory ganglion cells
 - amazing adaptation of brain \rightarrow
can analyze new sensory information by training
- Mimimum resolution:
 - project image divided into pixels onto normal retina
 - 25 x 25 \rightarrow can detect movement
 - 32 x 32 \rightarrow can read books with some help

BME2009

How many pixels are required?



8 x 8


32 x 32

320 x 320

- 100 pixel image (10 x 10)
- 625 pixel image (25 x 25) : enable mobility
- 1024 pixel image (32 x 32) : partially useful vision
- 10,000 electrodes (100 x 100) : ambitious goal

BME2009


Subretinal vs. Epiretinal



BME2009

Subretinal stimulator


- Advantage
 - use intraretinal neural network
 - Stimulation target : bipolar cell



BME2009


Optobionics Design

- **Used a simple photodiode array chip, no biasing, photovoltaic mode**
- Natural light is not strong enough to generate electrical power necessary.
- Experiments under FDA permission on 'in vivo biocompatibility'
- Operation on 12 patients with RP, from 45 to 76, in 2000.
- Nevertheless, some improvement of vision is reported in all patients.



Epiretinal stimulator


- Stimulation target: Retinal ganglion cell
- Doheny / NCSU / SecondSight



BME2009

Doheny / Second Sight


- 4 x 4 Pt electrode
- powered by Clarion® cochlear implant system
- Experiment results from patients since December, 2001
 - Detect light direction : 90-100%
 - Orientation of "L" : 75%
 - Location of white object on black background : 80%



BME2009

Optic nerve stimulator


- Stimulation target : Optic nerve
- Disadvantage
 - High density of axons (1.2million/2mm diameter) with tough dura
 - Uncertain topological arrangement of nerve fibers
- UCL (Universite Catholique de Louvain), Belgium
- First human trial in April, 1998
 - spiral cuff electrode (Au-Ti) with 4 electrodes
- Second human trial in 2003
 - 8 electrodes



BME2009

Visual cortex stimulator



- stimulation target : visual cortex
- Advantage
 - Broad applicability.
- Disadvantage
 - Complex topology of neurons.
 - Stimulation at a large area of visual cortex
 - Epileptogenic (간질발생)
- Utah, Dobelle Inst., Kresge



BME2009

Visual cortex stimulator

- Dobelle (1976):
 - **64 Pt electrode**, 8 x 8 array on 3 mm centers in Teflon ribbon cable matrix
 - use only 6 electrodes among them.
 - connection with the outside through penrose drainage
- Dobelle (2002):
 - Insisted that a operated patient drive and sense color



BME2009

Artificial retina : problems to be solved

- where should a stimulator be implanted?
- How can visual information be transmitted to the stimulator?
- How can energy be supplied to the stimulator?
- What is the optimal stimulation strategy?
- How can the eternal and best stimulator be made?


BME2009

SNU artificial retina

- Nano-bio system research center (for 9 years, 3 main subjects)
 - supported by KOSEF(KOREA SCIENCE AND ENGINEERING FOUNDATION)
 - development of electrode for retina stimulation as a subject in Neural chip/MEMS
- Nano artificial vision research center (for 6 years)
 - supported by Minister for Health, Welfare and Family Affairs
 - development of SNU artificial retina system and application to a human body
- Collaborators
 - ophthalmology
 - physiology
 - Electrical Engineering & Computer Science
 - biomedical engineering


BME2009

SNU artificial retina




BME2009

Retina Implant 人工網膜(視覺)



BME2009

Polyimide-based Retinal Prosthesis



Electrode spec.

Whole structure size : 3mm x 17.8mm x 16um
Size of each site : 50um x 50um
Number of sites : 25
Site spacing : 300 um
Insulation layer : Polyimide (PI2525)
Shape of head : Rectangular, Triangular
Head size: 3mm x 3mm



Acute stimulation
Experiment (rabbit)

BME2009

Micromachined Silicon Tacks

□ Fabrication


- material: silicon, oxide on silicon, parylene on silicon
- full scale dimensions: length 3 mm, width 300 μm, height 300 μm
- half scale dimensions; length 1.5 mm, width 150 μm, height 150 μm





BME2009

Surgical methods

□ subretinal insertion




□ epiretinal insertion




BME2009

Epi-retina implant using Polyimide electrode and Micromachined Tacks




BME2009

Setup for stimulation and recording



BME2009

Electrophysiological recording after insertion of retina stimulator




Visual evoked Cortical potential (VECP)

Electrically evoked cortical potential (EECP)

After optic nerve cutting


BME2009

Range of increase in glycometabolism of visual cortex (PET image)




Artificial Retina System for Animal Study

- Simultaneous Wireless Transmission of Data for Stimulation and Power (Class-E amplifier)




BME2009

Fabricated ICs



Data/Power Receiver Chip

- Rectifier/Amplifier
- Data Decoder (125Kbps)
- Voltage Regulator (5V)
- Internal Clock Generator(125KHz)




Current Stimulator Chip

- Current Source (7 channel)
- Biphasic Wave Form Generator

BME2009

Fabricated Artificial Retina System for Animal Study



BME2009